


EXPERIENCE THE COLOSSAL POWER & TORQUE OF THE SMALL AND MIGHTY ATOM



The Atom micro-head highspeed offers a lightweight and slim body design. Which allows for easy access to hard-to-reach areas with superior visibility and accessibility without compromising comfort. Delivering 16W of power the Atom allows you to cut smoothly and confidently without the loss of torque, elevating the Atom to the highest level of performance within it's class of micro-head highspeeds.



MICRO-HEAD

Micro head measuring 10.5mm by 10mm makes the Atom ideal for a smaller mouth.



LIGHTWEIGHT & SLIM BODY

Slim body design and weighing only 1.5oz the atom is designed to fit comfortably in all hand sizes.

ATOM

1 YEAR | ITEM# ATOM/ ATOM-FO

ATOM (NON-OPTIC)

\$399.99 | 1 Year Warranty | ITEM# ATOM

- Power - 16 Watts
- Ceramic Bearings
- Weight - 42.5g
- RPM - 430,000
- Head Height - 11.0mm
- Head Diameter - 10.2mm
- Fixed - 4 Hole

ATOM (FIBER OPTIC)

\$459.99 | 1 Year Warranty | ITEM# ATOM-FO

- Power - 16 Watts
- Ceramic Bearings
- Weight - 40.5g
- RPM - 430,000
- Head Height - 11.0mm
- Head Diameter - 10.2mm
- Fixed - 5 Hole Optic

THE MOST POWERFUL MICRO-HEAD AVAILABLE



ULTRA LIGHT, QUIET AND POWERFUL

Offering an impressive 16 watts of power allowing you to cut smoothly and confidently without any loss of torque



SMALLER MICRO PUSH BUTTON HEAD

Designed with a smaller micro push button head. Provides superior visibility and accesibility.